

## **RELAXATION THROUGH HEAT SPA**

Thanks to the selectable and soothing heat function, you can relax twice. With the carbon fibers equipped in the backrest, which serve as an infrared heater, the heat is distributed optimally and ensures relaxation of your tense muscles during the massage.

### INTELLIGENT 3D MASSAGE TECHNOLOGY

The 3D massage built into the Alphasonic III allows you to relax at the highest level. The complex and sophisticated massage system gives you the feeling of a real, human hands massage. Feel the depth and efficiency of the massage in the Alphasonic III.

# CALF KNEADING & FOOT REFLEXOLOGY MASSAGE

Calf kneading and foot reflexology massage are also part of the Alphasonic III's repertoire. Using compression massage, the calves are surrounded by airbags and pampered with air pressure. The foot massage of the Alphasonic III stands out from other models with the additional function "foot scraper". Here, the massage additionally affects the arch of the foot and the heel, so that after standing, walking or running for a long time, your feet can also relax completely.

# FEEL WEIGHTLESS WITH ZERO-G!

Using advanced Zero-Gravity technology, the Alphasonic III brings your entire body to one level, creating a feeling of weightlessness. In this position your heart, legs and back are on alignment. Due to the optimal weight distribution, you thus come to an immediate relaxation and an effective stress relief.




Smartphones that support the NFC function can easily be charged inductively with the charging surface of the AlphaSonic III.

WELLNESS COMFORT FOR YOUR HOME.



Dimensions upright	142 x 82 x 124 cm
Dimensions reclining	179-204 x 82 x 90 x cm
Weight	113 kg
Voltage	220-240V ~ 50/60Hz
Nominal power	150 W
Rated runtime	20 minutes
braintronics <sup>®</sup>	21 minutes
Certificates	@CEM@X



Width of seating area 53 cm | Width of shoulder area (inside) 58 cm

AVAILABLE COLOURS



Casada International GmbH

3104 Paderborn | Germany

Obermeiers Feld 3

ww.casada.com





V





# **SPECIALS**



# **BRAINTRONICS®**

Braintronics® technology synchronises your brainwaves through audio stimulation. Combined with a specially developed massage programme your stress-ridden body will be able to relax. No matter whether you'd like to relax, dream or learn – anything's possible.

ARTHROSE TRACTION

Arthrosis Traction consists of gentle stretch-

ing pressure. It has a stimulating effect on

the cartilages, which enhances the transport

of bodily fluids, nutrients absorbation and

If necessary, the double cushion can be

placed to reduce the neck and shoulder

can alleviate symptoms of arthrosis.



Wireless data transmission.

Key acupressure points are massaged.

# REFLEXOLOGY

-

Reflexology regenerates foot muscles and

### HUMAN HANDS FEELING

Specially designed massage heads imitate

With one touch the quick start buttons in your armrest allow you to control the most important functions.



functions

Wireless charging function for NFC-smartphones.

# MAGNETIC THERAPY FUNCTION

a unique massage function using magnet-

Kneading and pressure massage functions in the leg airbags ensure that any tension of the leg muscles that may be caused by long periods of standing are relieved


**BLUETOOTH** 

ACUPRESSURE POINTS

has a harmonizing effect on almost all or-gans inside your body.

the hand motion of a professional masseur The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra. QUICK START

### DUAL SENSOR BODY SCAN



Use the voice control to activate selected

WIRELESS CHARGER

Air bags for the left and right palm have ic field therapy.

CALF KNEADING



MaRM

0.0

00

### Pleasant alternation between relaxation and releasing the muscles - also perfectly suited for a final massage, for relaxation and recovery of the muscles.

# KNEADING

Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.

### **AIR COMPRESSION**

Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

# KNEAD & TAP (DUAL ACTION)

The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.

# TAPPING

Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.

### SOLE ROLLER

The foot massage stimulates the nerve endings and certain organs of the body, which also relaxes the muscles.

# L-SHAPE

Ô

You can expect a profound full body massage from the neck down to the buttocks.

# S-HNF

With this automatic program the massage heads are working very close to your body, they massage the whole back down to your buttocks. The S-shaped mechanic is constructed after the natural anatomy of the spine.

Ø

ZERO-G

**ĀUTO** 

MEMORY

For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!

# ZERO GRAVITY

ZERO SPACE

You put yourself in the position of an astronaut – that "weightless feeling".








Using the intelligent memory function the massage chair remembers your favourite massage exactly and begins immediately without a scan



An individual applicable to any size massage area.



The distance between the both massage



Further relaxation and fun with the music function

ADJUSTABLE FOOTREST The footrest can be extended, this way it is adaptable to different body sizes. ADJUSTABLE PILLOW





The adjustability of the backrest and footrest contributes to personal well-being of the body position.

 $\wedge$ 

 $\mathbf{i}$ 

3D

SHIATSU

# ADJUST SHOULDER

The massage heads can be moved upward or downward to a suitable position in the shoulder area.







Whole body airbag massage. Different areas are individually selectable.


### BACK STRETCH ROLLING Targeted stretch massages activate the body and have an invigorating effect on the





AIR ARFA





der/neck area – ideal for easing tension

# WAIST STRETCH

ly counteracts tensions.

The massage heads can massage a particu-















SHOULDER GRASP Targeted kneading massage in the shouland pain in this area.

Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effective-











The human brain is a complex organ consisting of millions of little nerve cells, called neurons. Each cell constantly produces tiny electrical impulses. Every thought, every sound and every impression, everything is transmitted and communicated by those small electrical impulses. Electricity is, so to say, the language of the brain. The sum of those impulses can be detected and diagrammed, e.g. by using an EEG. Charted on a graph this electrical activity forms a wave-like pattern, known as brain-

waves.

With braintronics®, combining audio stimulation and guided meditations with a soothing massage, these brainwaves can be stimulated in order to reach different states of relaxation in the brain. Whether you want to sustainably enhance your concentration, or just relax and pamper yourself - everything is possible with a massage chair from the Casada family.









App Store

- 🍉 Google play