HIGHLIGHTS



leasant alternation between relaxaion and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.



Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.



AIR COMPRESSION

Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



CARBON HEATING

Deeply penetrating Carbon infrared heat, blood vessels expand, the blood low increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

(NEAD & TAP (DUAL ACTION)

The techniques of kneading and tap-

oing are combined in order to relieve

rension and stimulate blood circulation



at the same time.

Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.



ve endings and certain organs of the body, which also relaxes the muscles.

You can expect a profound full body nassage from the neck down to the

The foot massage stimulates the ner-



With this automatic program the massage heads are working very close to your body, they massage the whole back down to your buttocks. The S-shaped mechanic is constructed after the natural anatomy of the spine.



For the lying-position the massage chair noves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!



You put yourself in the position of an astronaut – that "weightless feeling".



D MASSAGE

The intelligent 2D massage mechanism allows the massage chair to measure vour back to determine a more exact position of the sections of the spine.



AUTOMATIC

Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd



BACK AREA

An individual applicable to any size massage area.



ADJUSTABLE FOOTREST

The footrest can be extended, this way t is adaptable to different body sizes.



ADJUSTABLE PILLOW

If necessary, the double cushion can be placed to reduce the neck and shoulder massage intensity.



RECLINE POSITION

The adjustability of the backrest and footrest contributes to personal wellbeing of the body position.



ADJUST SHOULDER

The massage heads can be moved upward or downward to a suitable position in the shoulder area.



AIR INTENSITY Adjustable airbag massage intensity.

AIR ARFA Whole body airbag massage. Different areas are individually selectable.



Shiatsu (finger press) is a form of body herapy originating in Japan. The massage technique consists of soft, hythmic, far-reaching stretches and rotations. The aim of Shiatsu is to simulate individual parts of the body and to mobilise the muscles.





Nassage speed is adjustable.



The distance between the both massage heads can be set (narrow, medium



Further relaxation and fun with the music function



ULFTOOTH

wireless connection to your smartohone or tablet ensures a facilitated nusic experience.



ACUPRESSURE POINTS

Key acupressure points are massaged.



REFLEXOLOGY

Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.



HUMAN HANDS FEELING

specially designed massage heads mitate the hand motion of a professional masseur. The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical ands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra.



QUICK START

JSB-CHARGER

in your armrest allow you to control the most important functions.



While you relax, the integrated USB port can supply, for example, your phone with power.



CARBON LEG HEAT

s ideal for those who suffer from The heat dilates blood vessels and stimulates blood circulation. At the same time, muscle tone is relaxed, esulting in deep relaxation.





Seating size

Shoulder with

Nominal power

Rated runtime

44 cm

50 cm

78 ka

95 W


20 minutes


21 minutes

⊕C€EH[ℤ

220 - 240 V ~ 50/60 Hz

With one touch the quick start buttons





Obermeiers Feld 3 33104 Paderborn

Phone: 05254 / 930980 E-Mail: info@casada.com Web: www.casada.com



POLLUX

Maximum comfort for your home















EXPERIENCE THE ULTIMATE RELAXATION WITH THE POLLUX MASSAGE CHAIR -SPACE-SAVING AND COMPACT!

Don't miss the unique massage experience of the Enjoy the astronauts' technology and sink into com-Pollux massage chair. Space-saving with a comrange of different functions that will synchronise escape from your everyday life. your body and mind. Treat yourself with a choice The specially developed massage heads imitate of 6 automatic massage programs or select different massage techniques for an individual massage according to your needs. Due to its high quality How nice would it be to welcome your personal and long-lasting materials, the massage chair proves to be a constant and stylish mate. Benefit from its carbon lumbar heat for counteracting tense muscles.

plete relief for your back with the Zero-G function. pact design and at the same time offering a huge Feel weightless and gain the unique opportunity to

the feeling of a professional hand massage.

masseur whenever you want to in your cozy home?

